

# ScruffGuard

ScruffGuard enables you to handle and dose small animals gently, whilst avoiding bites and needlestick injuries.

ScruffGuard can be used when administering substances by subcutaneous, intraperitoneal and intramuscular routes. It can be used to implant microchips subcutaneously for identification purposes. A ScruffGuard is particularly helpful for trainees who are learning to handle small animals for the first time.



## Disposable Version

<b>Material</b>	Bagasse (Sugar Cane Pulp)
<b>Available Sizes</b>	Size 9
<b>Autoclave</b>	Max. 121°C
<b>Sterility</b>	Non-sterile
<b>Weight</b>	6g
<b>Colour</b>	Cream


## Reusable Version

<b>Material</b>	Nylon
<b>Available Sizes</b>	Size 6 - 10
<b>Autoclave</b>	Max. 121°C
<b>Sterility</b>	Non-sterile
<b>Weight</b>	20g
<b>Colour</b>	White

Product Code and Pack Size:

**20-SCRUFD25 ScruffGuard, disposable, pack 25**

**20-SCRUF9RE ScruffGuard, reusable, single pack**



**Website:** [www.datesand.com](http://www.datesand.com)  
**Telephone:** 0800 161 5831  
**International:** +44(0)161 274 1080  
**Email:** [sales@datesand.com](mailto:sales@datesand.com)



Available from:



[info@scanbur.com](mailto:info@scanbur.com)