



DOWNFLOW MODULE

IDEAL AS BUFFER STORAGE, STERILE
AND CAGE PREPARATION AREA



KEEP THE CLEAN EQUIPMENT CLEAN UNTIL YOU NEED IT

The downflow module allows you to create a sterile zone tailor-made to your requirements.


You decide how large you want the sterile zone to be. We join the modules on-site.

A quick, easy and cost efficient solution! Can you afford not to?



DOWNFLOW MODULE

IDEAL AS BUFFER STORAGE, STERILE
AND CAGE PREPARATION AREA



Ideal as:

- Buffer storage for autoclaved equipment (cages, lids, tops, bottles, etc.)
- Sterile corridor from the autoclave to e.g. preparation area, cage change area, etc.
- Cage preparation area for dispensing of bedding, diet, enrichment and water

The downflow module is a fully re-circulating unit providing HEPA filtered laminar air flow.

It is designed to protect equipment and procedures from contamination from the environment.

The downflow module operates at positive pressure.

For more information about **SCANBUR** products go to www.scanbur.com

info@scanbur.com