



DOWNFLOW

Module

Ideal as:

- Buffer storage
- Cage preparation area
- Sterile area



DOWNFLOW

Module

Keep the clean equipment clean until you need it.

The downflow module allows you to create a sterile zone tailor-made to your requirements. You decide how large you want the sterile zone to be. We join the modules on-site. A quick, easy and cost efficient solution! Can you afford not to?



Ideal as:

- Buffer storage for autoclaved equipment (cages, lids, tops, bottles, etc.)
- Cage preparation area for dispensing of bedding, diet, enrichment and water
- Sterile corridor from the autoclave to e.g. preparation area, cage change area, etc.

The downflow module is a fully re-circulating unit providing HEPA filtered laminar air flow.

It is designed to protect equipment and procedures from contamination from the environment.

The unit operates at positive pressure.

CONTACT US

Tel. +45 5686 5600
info@scanbur.com

SCANBUR 
Improving Life Sciences

www.scanbur.com

09.2014